

# Unlocking Soil Secrets: Mastering Precision and Bias in Metals, TCLP, Anions, and PFAS

**Craig Huff** Senior Technical Manager

#### Webinar survey – your feedback is appreciated!



### **Unlocking Soil Secrets**

- Key Learning Topics
  - Manufacturing processes and tips for optimizing digestion process
  - Soil method performance data
  - Anion testing
  - PFAS analysis methods
- Speaker Craig Huff
  - Senior Technical Manager

### Outline

- Metals in Soil Precision & Bias
  - Assigned values and homogeneity
  - Digestion hints
- TCLP Metals in Soil Precision & Bias
  - Assigned values and homogeneity
  - Observations
- Anions in Soil Precision & Bias
  - Assigned values and homogeneity
  - Extraction information & observations
- PFAS in Soil PT Studies Method Performance Observations





### Metals in Soil

- Data represent the last 10 years of PT studies
- Approximately 40 PT studies
- Represent thousands of data points per analyte
- Multiple types of soil utilized
- Represent multiple analytical methods
  - Data can be segregated by analytical methods, however it is consolidated in this presentation
- Represent multiple digestion methods
  - Unfortunately, the data could not be segregated by digestion method

#### Metals in Soil

**Precision & Bias Summary** 

#### Mean (% Recovery) Standard Deviation (%) Mean (% Recovery) Standard Deviation (%) Failure Rate (%) Failure Rate (%) **Historical** Analyte **Historical Historical** Analyte **Historical Historical Historical** Aluminum 86.4 24.1 3.2 Manganese 97.7 8.9 3.0 Antimony Mercury 13.7 5.0 43.4 57.7 1.6 98.2 Molybdenum Arsenic 86.2 9.2 2.2 83.4 10.6 2.4 Barium 94.4 8.8 2.7 Nickel 88.0 9.1 2.1 Beryllium 92.2 8.6 2.2 Potassium 13.5 93.8 3.0 Selenium 2.7 Boron 13.5 3.2 90.0 10.6 81.0 Cadmium 8.8 2.0 89.4 Silver 10.4 91.0 6.0 Calcium 95.2 Sodium 8.9 3.1 94.3 12.6 5.6 Chromium Strontium 9.5 90.2 9.3 2.1 94.9 2.9 Thallium Cobalt 90.5 8.6 2.2 88.5 9.9 2.6 Copper 90.6 8.5 2.3 11.0 4.1 Tin 89.4 1.6 Titanium 34.5 1.8 89.6 20.9 87.0 Iron Lead 9.3 2.3 Uranium 91.1 99.2 12.0 7.6 Lithium 89.1 19.0 25.0 Vanadium 86.6 11.2 1.9 Magnesium Zinc 9.7 93.2 11.7 2.1 93.8 2.1

Side note: AI & Fe have high native concentrations in the soil matrix

# Metals in Soil

# Waters<sup>™</sup> | ♦ ERA.

Assigned Values, Assurance of Homogeneity and Digestion Hints

- PT study assigned values = Study mean per TNI FoPT table
- CRM assigned values = ERA "made-to" values based on weights and measures + measured background concentrations
- Homogeneity testing is conducted by analyzing multiple samples, randomly selected from across the batch.
  - Samples are tested both across the batch and within each bottle before release
- Digestion Method = Method 3050B with ICP-AES and/or ICP/MS analyses
  - Mercury is analyzed per Method 7471(CVAA)
- Digestion Hints:
  - If using block digestors, make sure proper temperature is achieved, consistent and maintained in <u>each</u> well
  - Addition of HCI (per Method 3050), will enhance ICP-AES recoveries for many analytes!
  - Freshly prepared reagents are a must for mercury analyses!



### **TCLP Metals in Soil**

- Data represent the 5 years of PT studies
- 19 PT studies
- Multiple types of soil utilized
- Assigned values = ERA recoveries using Method 1311 extraction with ICP-AES and CVAA (Hg) analyses
- Failure rates based on acceptance limits equal to ± 3 standard deviations applied around the study mean for each analyte
- Only TCLP extraction fluid #1 was used

#### TCLP Metals in Soil Precision Summary

# Waters<sup>™</sup> | ♦ ERA.

	Average	Total	Average	
Analyte	Standard Deviation (%)	(n)	Failure Rate (%)	
Antimony	11.5	789	5.8	
Arsenic	11.2	1045	4.6	
Barium	10.0	1043	5.6	
Beryllium	11.3	756	4.9	
Cadmium	9.4	1067	3.3	
Chromium	11.0	1065	4.0	
Lead	16.1	1102	5.5	
Mercury	24.7	787	5.2	
Nickel	9.9	851	3.7	
Selenium	11.5	1037	3.4	
Silver	21.3	998	5.9	
Zinc	10.9	843	4.4	

#### ± 3 S.D. Acceptance Limits Also Note Hg & Ag %RSDs

# **TCLP Metals in Soil**

# Waters<sup>™</sup> | ♦ ERA.

Notes and Observations on ERA Internal and PT Data

- Percent recoveries using Method 1311 extraction yield only a fraction of actual amounts of each metal spiked onto soil
- Similar soil matrices (i.e., similar % sand/silt/clay content) yield consistent results based on ERA internally derived data and PT data
- Precision data should not be extrapolated to other real-world solid wastes or extractions using TCLP extraction fluid #2
- Extracting for shorter or longer periods than prescribed in Method 1311, will yield different recoveries (no surprise there...)
- Mercury and silver variability due inconsistent extraction efficiencies



### Anions in Soil

- Data represent the last 10 years of PT studies
- Approximately 40 PT studies
- Represent thousands of data points per analyte
- Multiple types of soil utilized
- Represent multiple analytical methods...but most common method is Ion Chromatography (Methods 9056x & 300.0)
  - Data can be segregated by analytical methods, however it is consolidated in this presentation
- Represents deionized water extraction only



#### Anions in Soil Precision and Bias Summary

Analyte	Mean (%) Recovery	Standard Deviation (%)	Failure Rate (%)	
	Historical	Historical	Historical	
Bromide	91.5	9.0	6.8	
Chloride	97.0	10.3	6.6	
Fluoride	33.6	26.5	3.3	
Nitrate as N	76.8	8.8	5.1	
Nitrate + Nitrite as N	76.8	8.8	3.2	
ortho-Phosphate as P	32.3	34.9	3.6	
Sulfate	85.2	14.0	7.1	

#### Note precision and accuracy for F<sup>-</sup> and PO4-P

### Anions in Soil



Assigned Values, Assurance of Homogeneity and Extraction Information / Observations

- PT study assigned values = Study mean per TNI FoPT table
- CRM assigned values = ERA "made-to" values based on weights and measures + measured background concentrations
- Nitrite is unstable in this matrix and does not contribute to Nitrate + Nitrite as N assigned value
- Homogeneity testing is conducted by analyzing multiple samples, randomly selected from across the batch
  - Samples are tested both across the batch and within each bottle before release
- Fluoride and Phosphate extraction efficiencies and overall precision can be improved by using weak HCI extraction solution...however, this solution cannot be used with this sample design due to method interference concerns and TNI FoPT table regression equations for precision

# **PFAS in Soil**

- Data represent the last 5 years of PT studies
- 20 PT studies
- Multiple types of soil utilized
- Failure Rates based on acceptance limits of 50-150%
- Most commonly reported methods:
  - DoD/DoE QSM table x
  - 1633 DRAFTx
  - 537(Mod)
  - ASTM D7968
  - EPA 8327
  - "Internal"



#### **PFAS in Soil** *Precision & Bias Summary*

Analyte	Mean (%) Study	Std Dev (%) Study	(n)	Failure Rate (%) Study
11chloroeicosafluoro3oxaundecane1sulfonic acid (11ClPF3OUdS)	91.7	18.3	308	3.0
9chlorohexadecafluoro3oxanonane1sulfonic acid (9CIPF3ONS)	96.2	13.5	304	1.0
4,8dioxa3Hperfluorononanoic acid (DONA)	89.4	16.4	313	1.0
Nethyl perfluorooctanesulfonamidoacetic acid (NEtFOSAA)**	87.5	17.5	93	0.9
1H, 1H, 2H, 2HPerfluorodecanesulfonic acid (8:2 FTS)	85.3	18.6	352	5.4
1H, 1H, 2H, 2HPerfluorohexanesulfonic acid (4:2 FTS)	88.4	14.4	331	1.2
1H, 1H, 2H, 2HPerfluorooctanesulfonic acid (6:2 FTS)	88.2	15.7	350	3.0
Hexafluoropropylene oxide dimer acid (HFPODA)	87.5	14.7	337	0.0
Nmethyl perfluorooctanesulfonamidoacetic acid (NMeFOSAA)**	85.7	15.2	91	2.0
Perfluorobutanesulfonic acid (PFBS)	91.5	12.4	384	1.8
Perfluorobutanoic acid (PFBA)	91.1	14.4	377	1.6
Perfluorodecane sulfonic acid (PFDS)	90.0	14.8	370	1.2
Perfluorodecanoic acid (PFDA)	90.2	13.7	377	2.5
Perfluorododecanoic acid (PFDoA)	91.5	13.7	376	1.9
Perfluoroheptane sulfonic acid (PFHpS)	91.6	12.6	347	2.3
Perfluoroheptanoic acid (PFHpA)	91.2	12.2	379	1.4
Perfluorohexanesulfonic acid (PFHxS)	89.5	11.0	380	2.1
Perfluorohexanoic acid (PFHxA)	90.8	13.3	380	0.3
Perfluorononane sulfonic acid (PFNS)	92.3	12.1	323	1.6
Perfluorononanoic acid (PFNA)	87.9	13.2	381	1.8
Perfluorooctane sulfonamide (PFOSAm)	92.5	16.2	351	0.5
Perfluorooctanesulfonic acid (PFOS)	93.1	14.1	389	2.4
Perfluorooctanoic acid (PFOA)	91.1	13.8	392	2.4
Perfluoropentanoic acid (PFPeA)	89.4	13.5	375	1.4
Perfluoropentane sulfonic acid (PFPeS)	89.2	13.6	331	1.2
Perfluorotetradecanoic acid (PFTDA)**	92.2	11.7	77	8.0
Perfluorotridecanoic acid (PFTrDA)**	92.6	13.0	120	3.8
Perfluoroundecanoic acid (PFUnDA)	92.8	13.9	376	0.7

#### \*\* Branched and Linear isomers. Also note low F.R's.

©2025 Waters Corporatior

# **PFAS in Soil**

# Waters<sup>™</sup> | �era.

Notes and Observations

- PT study assigned values = ERA "made-to" values based on weights and measures
- 4 Analytes composed of both branched and linear isomers
- Mean recoveries and standard deviations are actually pretty good overall
- PT Acceptance limits of ± 50% may be too wide?
- False positive reporting very insignificant! (data not illustrated in previous table)
- PFAS is a constantly evolving group of analytes. ERA's most recent studies contain an expanded analyte list of analytes and we anticipate this to change going forward



# Thank you!

# craig\_huff@waters.com





# Live Q&A Session: Unlocking Soil Secrets: Mastering Precision and Bias in Metals, TCLP, Anions, and PFAS

Please submit your questions into the Q&A Chatbox

https://www.eraqc.com